



## Reusable Products


**Perfusion Cannula**  
**Coronary Ostial Cannula**  
**Cardiac Suckers**

### Options

- Swivel or Fixed Tip Cannulas
- Full Range of Interchangeable Sucker Tips
- Handle and Tips Configured Per Your Application

### Benefits


- Non-Traumatic Blood Handling Capabilities
- Ergonomic Delivery
- Unique Leak Proof Swivel Tips
- Significant Cost Saving Over Disposables



Pioneer of Environmentally  
Friendly, Ergonomic and  
Economically Sound  
Surgical Devices

The unique design of Pemco perfusion cannulas allows user to combine different configurations of handles and tip styles with a luer or barb connection end. *View our full line of products at [www.pemcomed.com](http://www.pemcomed.com).*

- The unique ergonomic design of Pemco cannulas provides for proper placement of tip
- Pemco reusable cannulas are precision crafted and are available in Fixed or Leak Proof Swivel tip styles
- Tip styles and handle configurations shown are examples of full product line
- Manufactured from surgical stainless steel, Pemco cannulas can be sterilized using routine procedures.
- Long product life offers significant cost savings over disposables
- Refurbishing extends product life and provides additional cost savings



Handle Style:

P A A F T

## Handle Styles

(Can be configured with your choice of tip.)

1580-P-2-2MM Pediatric Handle

1580-F-7-5MM Flat Handle

1580-A-5-2MM Adult Handle

1580-T-8-4MM Tube Handle

1580-A-7-5MM Adult Handle



# Intracardiac Suckers Round Handle

Available in sets or individually

- Tips are grooved and slotted to allow for better flow and prevent clogging
- Suckers have 1/4 barbed end
- Alternate tip and handle styles are available
- Refurbishing extends product life and provides additional cost savings

## Round, Adult Handle Sets

The adult suckers are supplied as a set with four tips (#05, #06, #07 and #08). Handles and tips may be ordered as separate items. MODEL 1450: Intracardiac Sucker, Adult Handle at 30° Angle, Includes 4 tips. MODEL 1450-90: Intracardiac Sucker, Adult Handle at 90° Angle, Includes 4 tips.


### ADULT TIPS

#### Included with Set:

- 1450-05  
Intracardiac Sucker Tip  
10.5 mm diameter
- 1450-06  
Intracardiac Sucker Tip  
8.8 mm diameter
- 1450-07  
Intracardiac Sucker Tip  
8.0 mm diameter
- 1450-08  
Intracardiac Sucker Tip  
6.4 mm diameter

*Note: All above tips fit 1480-04 sucker handle.*

Adult Tips




1450 Intracardiac Sucker Set with four Tips



(Handle Only 1450-01)

1450-90 Intracardiac Sucker Set with four Tips



(Handle Only 1450-01-90)

10 1/2 inches long

## Round, Pediatric Handle Sets


The pediatric suckers are supplied as a set with two tips (#15 and #16). Handles and tips may be ordered as separate items. MODEL 1450JR Intracardiac Sucker, Pediatric Handle with 30° Angle, includes two tips. MODEL 1450JR-90 Intracardiac Sucker, Pediatric Handle with 90° Angle, includes two tips.

### PEDIATRIC TIPS


#### Included with Set:

- 1450-15  
Intracardiac Sucker Tip  
6.4 mm diameter
- 1450-16  
Intracardiac Sucker Tip  
4.0 mm diameter

Pediatric Tips




1450JR Intracardiac Sucker Set with two Tips



(Handle Only 1450JR-01)

1450JR-90 Intracardiac Sucker Set with two Tips



(Handle Only 1450JR-01-90)

7 inches long



# Intracardiac Suckers Flat Handle

Available in sets or individually


- Tips are grooved and slotted to allow for better flow and prevent clogging
- Suckers come with 1/4 barbed end
- Alternate tip and handle styles are available
- Refurbishing extends product life and provides additional cost savings

## Flat, Adult Handle

### MODEL 1480

Intracardiac Sucker  
Removable Tip, 7.0 mm diameter

(Note: All Adult Tips  
also fit 1480-04  
sucker handle)




1480 Adult Handle with One Removable Tip  
(1480-04 Handle Only / 1480-03 Tip Only)

## Flat, Adult Handle

### MODEL 1470

Intracardiac Sucker  
Fixed Tip, 4.7 mm diameter



1470 Adult Handle with Fixed tip

## Flat, Pediatric Handle

### MODEL 1460

Intracardiac Sucker  
Fixed Tip, 3.2 mm diameter



1460 Pediatric Handle with Fixed tip

**Pemco Reusable Instruments** reduce cost and save the environment by minimizing the use of disposables. Studies have demonstrated the following results for reusable instruments:


- Reduced disposable instrument purchases
- Reduced medical waste management expenses and scrapped plastics
- Reduced hazards to the environment, metals can be recycled
- Reduced risk to the health care worker, no need to “red-bag”
- Reduced licensing, handling, and inventory

### Other Benefits of Pemco Reusable Products:

- Proven safe and effective for over 50 years
- Routine instrument handling procedures
- Products are processed with your system
- Customized to your needs

**Pemco is committed to the best practices of your institution.**

3 UNITS PEMCO REUSABLE VS.  
200 CASES DISPOSABLE  
CANNULAS & SUCKERS



### Pemco Inc.

5663 Brecksville Rd. • Cleveland, Ohio 44131 U.S.A.

Phone: 216-524-2990 • Fax: 216-642-8646 • pemco@pemcomed.com • www.pemcomed.com

### European Authorized Representative

Medimark® Europe

Sarl 11 rue Emile Zola • 38033 Grenoble Cedex 2 • France